

Material
Processing
Technologies



Pharmaceutical Technologies International, Inc.
Princeton, New Jersey, USA

Particle Size Reduction MILLING TECHNOLOGY

Presented by:

Mina Ibrahim, P.Eng., MBA

Product Manager, Quadro Solids Division



Agenda



1. INTRODUCTION TO QUADRO, MPT, & IDEX
2. SIZE REDUCTION OVERVIEW
3. EVOLUTION OF MILLING TECHNOLOGY
4. OVERVIEW OF COMMON MILLING TECHNOLOGIES
5. CONICAL MILLING TECHNOLOGY
6. MILL SELECTION CRITERIA



About Quadro



- **History:** Since 1976
- **Manufacturing:** 45,000 ft² (4180 m²)
- **Employees:** 93
- **In-house Engineering**
- **Machining, Welding, Polishing, Electrical, Assembly**
- **ISO Registered, cGMP**
- **R&D Center**



About Quadro



- **Part of IDEX Corporation since 2007**
- **Member of Material Processing Technologies (MPT) platform along with**
 - **The Fitzpatrick Company, Illinois**
 - **Microfluidics, Mass.**
 - **Matcon, UK**



Size Reduction Overview



Advantages of Compressed Tablets (Oral Solid Dosage)

- Accurate dosage of medicament
- Easy to transport - bulk and by patient
- Uniform final product - weight and appearance
- Usually more stable than liquid medicines
- Release rate of drug can be varied
- Mass production - simple and quick & low cost



Size reduction is an essential process requirement in the practice of Solid Dosage Preparation

The capability to produce a tight particle distribution suitable for compaction and dissolution is directly dependent on the mechanism selected for size reduction

Size Reduction Overview



Tablet Manufacturing

- **Objectives**

- Uniformity
- Potency
- Batch to batch reproducibility
- Damage resistance
- Lack of defects

- **How**

- ✓ Powders must flow
- ✓ Powders must compress
- ✓ Particles must lock together
- ✓ PSD control → Weight control



Size Reduction Overview



Common Tableting Problems



Capping




Chipping



Sticking



Breaking



Discoloring



Porosity

Tablet weight is the key to controlling hardness and friability.

Controlling tablet weights within a tight range will contribute to better tablet hardness and friability.

Key weight control factors are product uniformity in particle size & density



Why Size Reduce

- Increase Surface Area
- Create Homogeneity
- Control Bulk Density
- Prepare Products for Post Processes

- Specifically for Tablets:
 - Increase bioavailability
 - Improve Flow
 - Reduce Segregation
 - Enhance Drying
 - Control Particle size
 - Repeatability – Batch to Batch

What Affects Size Reduction?

- Mechanical – Sizing Method (Type of Equipment)
- Fracture Mechanics of Particles – Types of Granules
- Properties of OSD ingredients:
 - **Active Pharmaceutical Ingredient (API)**
 - **Excipients - Inactive “helpers”:**
 - Anti-adherents/Lubricants: e.g. Magnesium Stearate
 - Binders
 - Wet: Gelatin, Starch, Sucrose, Glycol (dissolved in water or alcohol)
 - Dry: MCC, Polyethylene Glycol
 - Fillers: Lactose, Sorbitol, Calcium Carbonate
 - Flavouring/Colouring
 - Preservatives: Benzoic Acid, Cresol, Parabens, etc.
- Other physical properties – friability, toughness, abrasiveness, corrosiveness, etc.


Common Size Reduction Mechanisms

Size reduction equipment is available in many different designs, however, they all stem from four basic principles:


- **IMPACT:** particle concussion by a single force
- **COMPRESSION:** particle disintegration by two rigid forces
- **SHEAR:** produced by particle to particle interaction
- **ATTRITION:** arising from particles scraping against one another or against a rigid surface




EVOLUTION OF MILLING TECHNOLOGY



Stone Grinder




Roll Crusher



Lump Breaker



Hammer Mill



Oscillator



Tornado Mill



Conical Screen Mill

The most common method over
the last 30 years


Approx. 80 Years of
recorded history

Approx.
50 Years

Milling



One of the most essential process requirements in the practice of Solid Dosage Manufacturing




Common Milling Applications in the Manufacturing Process

- Dispensing – **De-agglomeration and security screen**
- Pre-Milling – **Particle Size Distribution**
- Post Granulating – **De-agglomeration/Dispersion**
- Dry Milling – **Sizing Dried Blend**
- Final Milling – **Size/De-lump/Calibrate**
- Reclaim - **Off-Spec Tablets/Compacts**

Milling



Typical Dispensing layout. Dust Free Design




Milling



Wet Dispersion Prior To Drying - Typical Integrated Design


High Shear Mixer/Granulator



Milling




Dry Milling After Fluid Bed Dryer Typical Integrated Design c/w Vacuum Transfer



Tablet Manufacturing



Wet Granulation





Comil U20

Wet Milling - Direct Discharge



**Fully Integrated Solid Dosage Preparation Plant
(Class 100,000 Room)**


A large industrial wet milling machine, the Comil U20, is shown in a clean, white laboratory or industrial setting. The machine consists of a tall, cylindrical stainless steel mill with multiple ports and a control panel. A flexible hose connects the mill to a discharge system. The discharge system includes a hopper, a grinding chamber, and a collection bin. The entire setup is mounted on a stainless steel frame with wheels. The background is a plain white wall with a ceiling light fixture.

Comil U20

Wet Milling - Direct Discharge

Hammer Mill




- High shear mechanism
- Various In-feed designs
- Variable speed, blade & hammer assembly
- 120° discharge area
- Common output range
(6" – 12" – 30" wide screens)



Hammer Mill Cont...

Blade & Screen Types

- The blade assembly is reversible.
- The most common blade arrangement is one blunt edge and one knife edge.
- Product can be hammered or cut.
- Hardened Blades available for abrasive applications



Hammer Mill Cont...

- Hammer Mills require control feed.
- Changes in feed rates may change product retention time.
 - will effect products that can easily dense
 - increased fines & friction
 - will effect products with low melting temperature



Hammer Mill Cont...

Advantages

- Wide range in Size
- Medium to High Shear
- Vertical/Horizontal Designs
- Blades/Screens Interchangeable
- Suitable for Milling Hard Materials

Disadvantages

- High Noise Levels
- % Fines High
- Must be control-fed
- Belt Slip Common
- High Dusting
- Ventilation Requirement
- Screen change complex
- Difficult to Scale-Up



Oscillator

- This machine was commonly used in the past for low shear applications.
- Some similar designs are continuous and do not oscillate.
- Suitable for low volume manufacturing.



Oscillator Cont...

- Uses mesh screens, not perforated plates.
- Cast Body
- Discharge - tray or drum
- High Wear rate.



Oscillator Cont...

Advantages

- Gentle
- Easy to operate
- Fixed speed
- Low cost equipment
- Low Tech Functions
- Portable

Disadvantages

- Low Capacity
- Metal to Metal contact
- Non GMP design
- Not suitable for integrated processes
- Cleaning - complex
- Loss of Active material


Conical Milling




- Infeed falls into conical screen chamber
- Rotating impeller imparts vortex flow pattern to infeed material
- Centrifugal acceleration forces particulates to screen surface
- Particles are continuously delivered to “action zone” between screen and impeller
- Particles are size reduced (as fine as 150 micron) and instantly discharged through screen openings



Conical Milling



**Underdriven Comil
(Invented 1990)**



**Overdriven Comil
(Invented 1976)**



Conical Milling – Quadro COMIL



Scale	Quadro Comil Model		Power	Standard Impeller speed	Scale-Up Factor	Tip Speed M/sec (Ft/min)	Screen Diameter	Capacity Lb/hr (kg/hr)
	OVERDRIVEN	UNDERDRIVEN						
Lab		U3	0.246 KW (0.33 hp)	4500 RPM	0.25X	14.2 (2800)	2.55" (65mm)	From 3oz/100g to 220lb (100kg)
		U5	0.375 KW (0.5 hp)	3450 RPM	0.5X	14.2 (2800)	3.25" (83mm)	425 (195)
Pilot	197	U10	1.5 KW (2.0 hp)	2400 RPM	1 X	14.2 (2800)	4.84" (123mm)	800-850 (360-390)
Production	194	U20	4.0 KW (5.4 hp)	1400 RPM	5 X	14.2 (2800)	8.2" (208mm)	3900-4250 (1750-1950)
	196	U30	7.5 KW (10 hp)	900 RPM	10 X	14.2 (2800)	12.17" (309mm)	7800-8500 (3500-3900)
Large Production	198		15 KW (20 hp)	450 RPM	20 X	14.2 (2800)	24" (609mm)	15,600 (7000)
	199		22 KW (30 hp)	360 RPM	40 X	14.2 (2800)	30" (761mm)	20,000 (9000)




Conical Milling – Quadro COMIL




OVERDRIVEN **QUADRO**® COMIL®-LAB TO PRODUCTION SCALABLE EQUIPMENT

Quadro Comils are completely customizable to meet your specific processing requirements



UNDERDRIVEN **QUADRO**® COMIL® –LAB TO PRODUCTION SCALABLE EQUIPMENT

Quadro Comils are completely customizable to meet your specific processing requirements




Conical Milling




Critical Factors for Optimum Conical Milling Characteristics

Close impeller / Screen Gap


- Less fines
- High Yield



Critical Milling Factors: Close Gap



Critical Milling Factors: Close Gap




Conical Milling



Critical Factors for Optimum Conical Milling Characteristics

Proper Tooling Selection – Screens


Effect of screen hole size on particle size distribution: generally a finer screen produces more fines and less overs.



Conical Milling




Critical Milling Factors: Proper Tooling Selection - Screens



R

R – Round holes
(Dry Material)




Q

Q - Square holes
(Wet Material)



S

S – Slotted
(Pseudo Plastic)



G


G - Grater holes
(Hard & Dry)



Conical Milling




Critical Milling Factors: Proper Tooling Selection - Impellers



1=Highest 4=Lowest




Comil Impellers – Round Arms #1601




Round arms - primarily for dry sizing, some wet milling

Comil Impellers – Rectangular Arms Positive Leading Edge #1607



Square Arms – “Universal” for wet milling and dry sizing

Critical Milling Factors: Screens - Apparent Hole Size



- Accepted definition of Fine Milling is **psd** between 5 - 100 Microns and for Micronization **psd** between 1 - 30 micron in diameter.
- It is possible to use some of the previously discussed equipment to reduce the particle size distribution of a product down to this range (Hammer Mill) however, distribution curve can be fairly wide spread and possibly even bimodal whereas a tight **psd** and unimodal curve is the goal of most processes.
- Equipment commonly used for fine milling are: Pin Mills, Hammer Mills, Fine Grind, & Jet Mills

Fine Milling



SIZE REDUCTION CAPABILITY COMPARISON


Comil																	
F10 Fine Grind																	
Hammermill																	
Pin Mill																	
Jet Mill																	
Micron	-5	-2.5	1	5	10	25	38	45	75	125	150	180	250	300	425	600	1000
US Mesh	-	-	-	-	-	-	400	325	200	120	100	80	60	50	40	30	18



Fine Milling



Quadro Fine Grind F10



1. Collector Cover
2. Product Collector Body
3. Product Hopper
4. Pneumatic Vibrator
5. Outlet container
6. Rotation Hinge
7. Milling Head
8. Screw Feeder
9. Control Panel
10. Access Panel



Fine Milling



- Fine Grind F10 was developed to produce tailored **PSD** between 15 and 100 microns.
- Mobile, stand alone system (a complete plant) operates at low noise, dust heat and energy consumption.
- **The operating principle;**
 - control feed product into upper conical screen chamber.
 - a rotating impeller calibrates incoming material.
 - calibrated product then passes through to the lower chamber
 - a second intensifying impeller accelerates the particles.


Upper Chamber



Lower Chamber




Fine Milling



F10 Breakthroughs


- Very tight Particle Size Distribution
- Very high Product Recovery rate (>99%)
- Dust-tight
- Mobile, All-In-One unit, no ancillaries required
- Operator-friendly: Easy to clean & Low-Noise Operation
- Sanitary, GMP design; Developed specifically for Pharmaceutical API Industry

Technology Comparison – Lactose 200M



Technology Comparison – Noise

Sound Level (No load / 1m away)



Case Study – Teva, Israel

MILLING EXPERIMENTS WITH PROTOTYPES OF THE QUADRO FINE-GRIND F-10 MILL

ABSTRACT

Quadro Fine grind (F-10) is a versatile mill which is effective in the range of 20-60 μm . Three Prototypes were received for experiments in Teva's API division. Different active pharmaceutical ingredients were tested to explore the mill's performance. Based on the conclusions, the final model was built. Improvements included enlargement of the milling chamber, Control over speed and vacuum, and introduction of water cooling. It was found that the milling range covers the particle size reduction range obtained today by either single or multiple milling in hammer-mills, and can provide comparable results to those of a pin-mill.

Case Study – Teva, Israel

INTRODUCTION

Teva's API division manufactures over 200 molecules for various pharmaceutical clients. The physical properties of the products are tailor-made in order to meet various customer requests and optimize the formulation [1]. Because of the large number of products and different physical grades, it is required that mills will be versatile, i.e. capable to produce a wide spectrum of P.S.D by changing only the operating parameter.

The P.S.D range of ~20-40 microns is considered to be difficult to obtain. Larger particles can be controlled by Hammer-mills, Comils or other mechanical mills. Particles under 20 microns can be obtained by fluid-jet mills. However, only few mills can obtain narrow P.S.D in this range without having too many fines or oversized particles. One of these mills is the Pinmill [2]. Few main drawbacks of this mill are the heat generation and the very narrow gap that make it prone to blockages. Therefore, a great interest was found in the Quadro Fine grind (F-10). Two prototypes were tested, and based on Teva's findings, the final version was constructed and successfully applied in routine production.

Case Study – Teva, Israel

Teva Paper at CHoPS Conference Italy, Aug 2006
“Development of the F10 in Teva, API”

Paper Synopsis

Goal: PSD 20 to 40 μm range

Previous: Pin Mill. Heat changed product characteristics.
Narrow gap between pins prone to blockage
(9 hrs to clean vs. F10 at 1 hour)

Validation: 6+ API's validated with F10

Case Study – Teva, Israel

Customer Requirement	Observations & Discussion:			
F10 comparison versus Pin Mill and Hammermill	Material “A”: F10 vs. Pin Mill / Hammermill			
	Material A	D ₁₀ (µm)	D ₅₀ (µm)	D ₉₀ (µm)
	Unmilled	60	180	410
	Pin Mill	2	15	45
	Hammermill Double Pass	4	20	50
	F10 Single Pass	1.6	11.9	49.4

1. PSD Comparison between F10 versus Pin Mill and/or Hammermill results provided equal or better PSD distribution.
2. Material “A” is a proprietary pharmaceutical API

Case Study – Teva, Israel

Customer Requirement	Observations & Discussion:			
F10 comparison versus Hammermill	Material “B”: F10 vs. Hammermill			
	Material B	D ₁₀ µm	D ₅₀ µm	D ₉₀ µm
	Unmilled	12.73	66.33	211.83
	Hammermill	8	50	150
	F10	3.44	18.69	63.33

1. PSD Comparison between F10 versus Hammermill results provided better PSD distribution.
2. Comil was also tested: D90 180 µm, D50 70 µm, D10 10 µm
3. Material “B” is a proprietary pharmaceutical API

Case Study – Teva, Israel

Customer Requirement	Observations & Discussion:			
F10 comparison versus Hammermill	Material “C”: F10 vs. Hammermill			
	Material C	D ₁₀ µm	D ₅₀ µm	D ₉₀ µm
	Unmilled	24.33	118.91	339.14
	Hammermill	7.96	57.34	157.62
	F10	7.59	30.84	85.04

1. PSD Comparison between F10 versus Hammermill results provided equal or better PSD distribution.
2. Material “C” is a proprietary pharmaceutical API

Case Study – Apotex, Canada

Customer using Hammermill: 4-5 passes for $d_{90} = 70 \mu\text{m}$

F10: $d_{90} = 53.6 \mu\text{m}$ (single pass) 7200RPM and $20.4 \mu\text{m}$ 8400RPM

Alendronate Sodium Trihydrate		Impeller Speed = 7200 rpm		Impeller Speed = 8400 rpm	
		Run 1	Run 1.1	Run 2	Run 2.1
PSD	Starting Material	PSD Run 1	PSD Run 1.1	PSD Run 2	PSD Run 2.1
D(v,0.1)	8.847 μm	3.503 μm	2.523 μm	2.694 μm	2.876 μm
D(v,0.5)	49.214 μm	18.03 μm	7.408 μm	7.585 μm	7.05 μm
D(v,0.9)	262.787 μm	53.601 μm	19.442 μm	20.451 μm	14.805 μm

First Pass

Second Pass

Typical F10 PSD Graph – MCC

Specific Surface Area:
0.275 m²/g


Surface Weighted Mean D[3,2]:
21.805 μ m

Vol. Weighted Mean D[4,3]:
42.411 μ m

d(0.1): 10.560 μ m




d(0.5): 31.714 μ m

d(0.9): 89.907 μ m



**F10 was run at standard speed (7200RPM), 045R screen;
Malvern Mastersizer 2000 Results**

Jet Mills & Micronizers



- The principle of micronizing fluid energy mills (also known as jet mills or spiral mills) is the size reduction of particles through inter-particulate collisions combined with surface collisions due to acceleration of product.
- These mills use accelerated fluid streams (normally compressed air, super heated steam or inert gas) to generate a high speed vortex which the particles are introduced into.
- The vacuum created by a venturi-nozzle propels the product throughout the milling chamber, forcing particles to collide with themselves as well as the chamber walls.

Key Components and attributes that affect micronization:

- **Nozzle design and direction of air jets**
- **Efficiency of air compressors**
- **Efficiency of filters and separators**

Mill Selection Criteria



- Properties of Feed Material:
- Size
 - Shape
 - Moisture content
 - Physical and chemical properties
 - Temperature sensitivity
 - Grindability
- Final Product Specification:
- Size
 - Particle size distribution
 - Shape
- Versatility of Operation:
- Change of speed and screens
 - Safety features



Mill Selection Criteria



Scale-Up:

- Capacity of the mill
- Production rate requirements

Dust Control:

- Loss of costly drugs
- Health hazards
- Contamination of plant
- Safety

Sanitation:

- Ease of cleaning and sterilization
- Design and material finish

Auxiliary Equipment:

- Cooling system
- Dust collectors
- Forced feeding



Mill Selection Criteria



- Economical Factors:
- Equipment cost
 - Power consumption
 - Space occupied
 - Labor cost



Ability to handle dust explosions

General guidelines for inert milling:

Minimum Ignition Energy: (ref. BS5958 Part 1; 1991)

< 500 mJ Low sensitivity to ignition. Solution: **Earth plant.**

< 100 mJ Recommended at this point that customer seek expert advice. Common solution: **Earth personnel.**

< 25 mJ Majority of incidents occur when MIE is at or below this level. Solution: **Inert with nitrogen.**

< 10 mJ High sensitivity to ignition. Solution: **Inert with nitrogen and monitor allowable oxygen levels.**



Thank you

